Hydrogenations with the new bpc 2 precise, reproducible, reliable, safe

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- unique volumetric gas dosing system
- proven technology
- turnkey solution – plug and play
- multiple built in safety features
- easy operation

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Facts and figures

- dosing of (hydrogen) gas at 1 – 140 bar
- control mode:
 - constant reaction pressure
 - constant flow rate
- continuous measuring of H₂ uptake
- gas consumption in mol, ml
- accuracy > 99%
- flowrates 0 – 225 liters n / min
- functional for:
 - H₂ (Hydrogen)
 - C₂H₄ (Ethylene)
 - CO₂ (Carbon dioxide)
 - CO (Carbon monoxide)
 - C₃H₆ (Propylene – soon)
 - O₂ (Oxygen – soon)

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The new bpc 2 and Büchi pressure reactors – the ultimate solution for hydrogenation

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Built on the experience and the outstanding performance of the first generation hydrogenation system, the bpc 2 offers unmatched reproducibility, proven accuracy, an intuitive interface, and operates at the highest level of safety.

Features

- reactor (reaction) pressure kept stable during entire experiment
- extremely precise with accurate reproducibility using our unique volumetric dosing
- continuous logging of consumed gas (to PC , USB memory stick)
- touch-screen operation, data display
- various internal safety interlocks
- automatic switch-off in event of leakage
- no need for on-site calibration
- sequence program, automation

Easy operation

1. purging with inert / active gas
2. set regulated pressure > delivery pressure
3. set delivery pressure and start

Simple installation

- 2 lines-in (Gas-supply)
- 1 line-out to reactor
- 1 vent line

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
Setup

bpc 2 controller

- touch-screen operation
- USB data memory port
- RS232 port
(for optional bls software)
- Ethernet

bpc 2 dosing unit

- volumetric gas dosing unit (VD*)
- pressure control devices
- safety interlock valves

Options

- measure / display of temperature, pressure, stirrer speed, alarm functions
- Sequence program (venting valve included): Example for system preparation before run:
 1. N2 purging
 2. Leak Test
 3. H2 purging
- bls software for bpc 2 remote control and data storage from and to PC

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
Easy and safe operation – valuable information

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The new bpc 2 is suitable for Büchi pressure reactors of different volumes, pressure ranges and materials.

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Applications

- Hydrogenations (catalytic)
- chemical research
- process development
- scale up

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Precise over wide ranges of operating pressure and flow rate

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max. supply pressure: 200 bar / 最大供給圧力:200bar

* delivery pressure is equal to reaction pressure / 吐出圧力は、反応圧力

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Please contact us for the 400 bar delivery pressure version

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
 **büchiglasuster**®
switzerland Pilot Plant and Reactor Systems